

# AbZelectPRO™

## Achieving high-titer productivity & accelerated timelines with enhanced cell line platforms




Brett Verstak\*, Blanca San Miguel, Oscar Swindley, Nathan Courtier, Kean Huedepohl, Ned Swift, Miles Wellington, Sophie Douglas, Noah Collins, Elizabeth Willows, Tania Fisher, Kalpana Wood, Tanya Knight, Beata Blaszczyk, Nicole Wakes, Rob Holgate  
Abzena, Babraham Research Campus, Cambridge, UK | \*email: brett.verstak@abzena.com

### Abstract

Our fully integrated AbZelectPRO™ platform enables the rapid delivery of stable, high-producing cell lines to support the production of complex biologics. The platform combines our AbZelectPRO™ CHO-K1 and AbZelectPRO™-KO GS knockout cell lines together with ProteoNic's 2G UNIC® premium vector technology and an optimised process to generate fast doubling, higher-producing, stable cell lines expressing up to 10 g/L of product, and DNA to RCB in as little as 10 weeks.

The AbZelectPRO™ platform supports a diverse range of modalities from antibodies and more difficult-to-express proteins such as fusion proteins, bispecifics, vaccines and other novel modalities. Furthermore, the platform also allows for program derisking by using early material generation from stable pools to remove activities such as downstream method development and formulation development from the critical path. Backed by our comprehensive analytics portfolio, the AbZelectPRO™ platform simplifies the development process and helps to deliver therapeutic proteins successfully to the clinic.

### Flexible options




**Figure 1. Multiple Cell Lines – One Platform.** Abzena integrates leading partnered technologies such as 2G UNIC® (ProteoNic) and CHOSOURCE™ GS KO cell line (Revvity) to streamline and accelerate CLD workflows. This unified approach enhances productivity, flexibility, and performance across platforms.

### Optimised CLD process



- High productivity** from stable pools and clones
- Optimized platform** that can accommodate **diverse formats**
- Fast doubling times** - DNA to RCB in 10 weeks
- High clonal stability** over 60+ generations
- Reproducible PQA** at all stages


The **CLD process has been optimised** to maximise efficiency, ensuring high productivity from stable pools and clones while providing a **robust process from DNA to RCB and beyond**.



**Figure 2.** Stable pool material generation (≥ 10 g for standard mAbs). Purified using state-of-the-art ÄKTA™ systems yielding ≥ 95% purity and low endotoxin.

**Figure 3.** Cyto-Mine® provides early population enrichment and single cell encapsulation for initial evidence of monoclonality with secondary evidence provided using Cell Metric® images throughout the single cell outgrowth time-course

### State-of-the-art analytical support



### Overview

**Trusted technology:** AbZelectPRO™ is built on an established CHO-K1 platform – well recognised in industry and by regulators.

**Experienced team with proven track record:** CLD team has delivered 150+ CLD programs with CHO-K1 cell line.

**Compliance assured full traceability:** Established robust processes compliant with ICH-Q5 guidance.

**Integrated platform** with end-to-end support under one organisation across all stages of the drug development process.